

FIXED / CERAMIC CROWN & BRIDGE



PREMIER DENTAL LABORATORY, L.C.
 "Committed to the Enhancement of Your Dental Practice"
 955 Dairy Ashford, Suite 222 • Houston, Texas 77079
(281) 589-6793

FOR OFFICE USE ONLY

CASE NO.	
LOG DATE	
MOD	
TRM	
WAX	
FIN	
POR	
QC	

PATIENT: _____ AGE _____ M / F

DOCTOR: _____

PREP _____ Tryin
 DATE _____ Finish DUE: **5pm**

Allow 10 days in lab plus 2 pu/del days (items marked * require 13 days in lab plus 2 pu/del days)
 Call Lab for due date on 8+ unit cases and custom shade cases.

FINAL SHADE: _____ STUMP SHADE: _____

- OCCL. STAIN:**
 None
 Light
 Medium
- INCSL. TRANSL:**
 Less
 Per Tab
 More



IMPLANT CROWNS:

- Screw-Retained
- Stock Abutment
- Modify Abutment
- Impression of Abutment
- "UCLA" type Custom Castable
- Telescopic Substructure
- Atlantis Abutment (Titanium / Zirconia)
- Procera Abutment (Titanium / Alumina / Zirconia)

Implant Parts Enclosed: _____

PORCELAIN to METAL:

- Porcelain fused to Non-Precious (Economy)
- Porcelain fused to Semi-Precious (Noble White)
- Porcelain fused to Precious (High Noble White)
- Porcelain fused to Precious (High Noble Yellow)
- Captek Alloy * / Bright Gold Alloy

VENEERS:

- LAYERED EMPRESS Esthetic Veneer*
- THIN VENEERS 3/10ths (e.max 400 MPa)

ALL PORCELAIN CROWNS:

- EMPRESS / e.max Pressed & Layered
- e.max (400 MPa) Pressed Full Contour

MARGINS:

- Metal Lingual Band
- Metal Collar (360°)
- Full coverage (no band)
- Porcelain butt shoulder

OCCLUSIONS:

- Porcelain Occlusion
- Metal Occlusion
- Metal Island (if needed)
- Full Metal Lingual
- Adjust Opposing

ZIRCONIA CROWN / BRIDGE:

- Layered Zirconia (CERCON / Lava)
- Full Countour ZIRCONIA HT (Cercon HT/BioZir)

FULL CAST:

- 50% Gold Crown / Inlay / Onlay (Noble Yellow)
- Economy Gold Crown / Inlay / Onlay (Noble Yellow)
- Full Cast Crown Non-Precious
- Full Cast Semi-Precious (Noble White)

INLAYS / ONLAYS:

- Porcelain Inlay / Onlay (e.max 400 MPa)
- Nano-Hybrid Premise Direct

PONTIC DESIGNS:



RADICA PROVISIONALS

DIAGNOSTIC WAX UP

STENT FOR PROVISIONALS

Special Instructions: _____
 (Tooth numbers)

PLEASE EVALUATE MY PREPS AND IMPRESSIONS

CALL ME

LIC No _____ SIGN. _____ D.D.S.